

3M™ Cavilon™ Durable Barrier Cream Application Guide



1. Squeeze a small amount of Cavilon Durable Barrier Cream into the palm or on the fingertips of your gloved hand. (Start with a circle of cream about equal to the size of a quarter.)

More cream is not better.
Avoid using too much!

If using the sachet packet, squeeze the entire contents into the palm or on the fingertips of your gloved hand.



2. Apply cream to all skin areas exposed to urine or stool. Gently massage cream into skin until it disappears.
3. Reapply every 24 hours if there are 3 or fewer incontinent episodes per 24 hours.
4. Reapply every shift if there are more than 3 incontinent episodes per 24 hours. Suggested protocol: apply at the first incontinent episode of the shift.

Use 1 application per shift



Use 3M™ Cavilon™ Durable Barrier Cream for:

- To help with routine Incontinence-Associated Dermatitis (IAD) prevention for intact skin (urinary incontinence and/or occasional fecal incontinence episodes)
- Moisturization

Use 3M™ Cavilon™ No Sting Barrier Film for:

- To help prevent IAD and to protect skin when there are frequent episodes of mixed incontinence or frequent/continuous diarrheal stools
- Can be used on intact or damaged skin
- Moisture and friction protection (pressure ulcer prevention program)
- Adhesive trauma protection
- Protection from friction
- Peristomal skin protection

Important Note: Refer to the product packaging for complete information related to the use of Cavilon Durable Barrier Cream. This application guide is intended only as a quick reference guide.



Critical & Chronic Care Solutions Division

3M Health Care

2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-228-3957
www.3M.com/Cavilon

3M and Cavilon are trademarks of 3M.
Please recycle. Printed in U.S.A.
© 3M 2013. All rights reserved.
70-2011-5472-4