

Easy to Start, Easy to Adjust for Optimised Outcome

Easy to Start

- 01 Calculate initial settings based on guidelines¹ recommendations using pump total daily insulin (TDI), if transitioning from multiple daily injections, calculate pump TDI by using pre-pump TDI x 0.75 or use Profile Settings Calculator (when available).

Basal Rate

50% Pump TDI / 24 hours

Correction Factor

95 / Pump TDI (mmol/L)
1700 / Pump TDI (mg/dL)

Carb Ratio

450 / Pump TDI


- 02 Enter weight and pump TDI to turn on automation

Easy to Review


✓ Assess glycaemic goal achievement

- + Review glycaemic goals and assess if personal glycaemic goals are being met by reviewing the time in range comparison
- + NOTE: CGM-based targets are general guidelines. Individual needs vary.²

CGM-Based Targets





Time in range comparison




✓ Evaluate time in automation, insulin delivery, and bolus behavior

- A Time in Automation:** Goal >90%
- B Sleep Activity:** Goal 7 times a week
- C Average Daily Dose:** Use if settings require recalculation
- D Bolus Behavior:** Assess manual bolus vs. AutoBolus*



✓ Identify patterns of hypoglycaemia and hyperglycaemia



Easy to Adjust

- ✓ Optimise settings when treatment goals aren't met
- ✓ Adjustable settings: Basal rate, Correction Factor, carb ratio
- ✓ When in doubt, prioritise **Correction Factor** to optimise both basal modulation and correction boluses

Stronger Correction Factor to maximise TIR without significantly increasing TBR^{3†}
 < 88 / pump TDI (mmol/L) or < 1580 / pump TDI (mg/dL)

Basal

Basal Modulation

Frees constraints and allows more dynamic basal delivery

Bolus

Strengthens Bolus

AutoBolus
Manual Correction Bolus

Smaller number = Stronger corrections to lower glucose levels

✓ Other Adjustments to Consider

Hypoglycaemia or hyperglycaemia not related to bolusing?



Change **basal rates**

Not bolusing consistently?



Strengthen **Correction Factor**

Hypoglycaemia or hyperglycaemia after bolusing?



Adjust **Correction Factor** and/or **carb ratio**

* If glucose values are predicted to be above 10.0 mmol/L (180 mg/dL), Control-IQ+ technology calculates an AutoBolus using the Personal Profile settings and a target of 6.1 mmol/L (110 mg/dL) and delivers 60% of that value. † With optimised Correction Factor based on "rule of less than 1600".

References: 1. Grunberger G, et al. *Endocr Pract.* 2014, 20(5):463-89. 2. Battelino T, et al. *Diabetes Care.* 2019;42(8):1593-1603. doi: 10.2337/dci19-0028. 3. Messer LH, Breton M. *Diabetes Technol Ther.* 2023;25(12):877-882.



Safety info: tandemdiabetes.com/safetyinfo.

© 2026 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, Tandem logos, Control-IQ, and Control-IQ+ are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. All third-party marks are the property of their respective owners. MAT-0664_2



TANDEM
Diabetes Care

tandemdiabetes.com