


NHS Forth Valley



HOW TO INSERT YOUR SLIM TUBE HEARING AID

Patient Information
Audiology Department



NOTES

How to insert your slim tube hearing aid

1. Place the hearing aid behind your ear



2. Hold the bend of the tube between your thumb and index finger with the dome pointing towards your ear canal.



3. Gently push the dome into your ear until the tube is flush against your ear.



4. Bend the sports lock back towards the rear of your ear until it curls neatly inside.



This is important to keep the aid securely in your ear.

Hearing aid inserted correctly:



5. To remove the hearing aid gently pull the bend of the tube until the tube is removed completely.



If you experience any problems, please DO NOT attend the department. Our open access clinic is CLOSED due to Covid 19 restrictions.

Contact can be made by phone or e-mail:

01786 434171

fv.audiology@nhs.scot

NOTES

We are happy to consider requests for this publication in other languages or formats such as large print.

Please call **01324 590886 (9- 5)**

or email fv.disabilitydepartment@nhs.scot

For all the latest health news visit www.nhsforthvalley.com

follow us on  [@NHSForthValley](https://twitter.com/NHSForthValley)

or like us on  at www.facebook.com/nhsforthvalley

If you can't go let us know!

Every month around 2,000 people across Forth Valley fail to turn up for hospital appointments. This costs the NHS millions of pounds each year and increases waiting times.

So if you are unable to attend, or no longer require your hospital appointment, please let us know so we can offer it to someone else.

SMOKING IS NOT PERMITTED ON NHS FORTH VALLEY PREMISES

This includes corridors, doorways, car parks & any of our grounds. If you do smoke on NHS premises you may be liable to prosecution and a fine.

NHS Forth Valley Administration Offices
Westburn Avenue, Falkirk, FK1 5SU
www.nhsforthvalley.com

