

How do I upload my Medtronic pump?

Before any appointment you should upload your data from your pump at home

Medtronic

To upload your data from your pump you will need to connect the USB part of your Contour next Linkmeter to your computer and have your pump close by. This link will explain how to do this:

[CareLink Personal Uploader | Medtronic Diabetes \(medtronic-diabetes.co.uk\)](https://www.medtronic-diabetes.co.uk/carelink-personal-uploader)

You will need to set up a carelink account if you don't have one. Here is the link:

[CareLink™ Personal Software – For People With Diabetes \(minimed.eu\)](https://www.minimed.eu/carelink-personal-software)