

How do I upload my Tandem T:slim X pump?

Before any appointment you should upload your data from your pump at home.

To download your data from your T:slim pump you are required to upload your data from your pump onto Glooko .

Step 1: Set up an account Link to our practice ID to share your data with us.

The clinic ID is 20-63945.

Step 2: Once you have a Glooko account and are logged in download your data from your Pump to Glooko by connecting your T slim to your computer/laptop/tablet with the USB cable.

*If you have an ipad you can buy a lightning or USB-C (ipad charge point) to micro USB adaptor for less than £10